X-RAY VIEW

Surface mounted luminaire for X-ray viewing. A lumieo opal diffuser material with 73% transmission and excellent lamp obscuration is installed into an aluminium frame.

The X-ray panel is wired with 3-lamp electronic ballast. The hotspot, if required, is obtained with a 15W CFL which is coupled to a push button switch.

665mm | 3x18w | 4x18w - T8 1270mm | 2x36w | 3x36w | 4x36w - T8 1570mm | 2x58w | 3x58w | 4x58w - T8

Applications: Hospitals, Doctors Rooms

MED-XR







